

UBE2U Polyclonal Antibody

Catalog No. E-AB-62576

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human UBE2U (NP_689702.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IHC 1:50-1:100

Data



Immunohistochemistry of paraffin-embedded Human liver using UBE2U Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse brain using UBE2U Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

UBE2U (Ubiquitin Conjugating Enzyme E2 U) is a Protein Coding gene. Diseases associated with UBE2U include Brachydactyly, Type E2. Among its related pathways are Remodeling of Adherens Junctions and Innate Immune System. Gene Ontology (GO) annotations related to this gene include ligase activity and ubiquitin protein ligase activity. An important paralog of this gene is UBE2B.

For Research Use Only