

VEPH1 Polyclonal Antibody

Catalog No. E-AB-63175

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human VEPH1 (NP_078897.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using VEPH1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:105kDa

Calculated Mw:19kDa/23kDa/89kDa/94kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

VEPH1 (Ventricular Zone Expressed PH Domain Containing 1) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include binding.

For Research Use Only