

GAS6 Polyclonal Antibody

Catalog No. E-AB-62753

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human GAS6 (NP_000811.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:50-1:100

Data




Western blot analysis of extracts of various cell lines using GAS6 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:69kDa

Calculated

Mw:41kDa/44kDa/68kDa/74kDa/79kDa



Immunohistochemistry of paraffin-embedded Mouse lung using GAS6 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene encodes a gamma-carboxyglutamic acid (Gla)-containing protein thought to be involved in the stimulation of cell proliferation. This gene is frequently overexpressed in many cancers and has been implicated as an adverse prognostic marker. Elevated protein levels are additionally associated with a variety of disease states, including venous thromboembolic disease, systemic lupus erythematosus, chronic renal failure, and preeclampsia.

For Research Use Only