

MEMO1 Polyclonal Antibody

Catalog No. E-AB-63109

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human MEMO1 (NP_057039.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IHC 1:50-1:100

Data



Immunohistochemistry of paraffin-embedded Human prostate using MEMO1 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Rat brain using MEMO1 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

MEMO1 (Mediator Of Cell Motility 1) is a Protein Coding gene. Diseases associated with MEMO1 include Periampullary Adenoma. Among its related pathways are Signaling by GPCR and Signaling by ERBB2.

For Research Use Only