

ARID2 Polyclonal Antibody

Catalog No. E-AB-62580

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant protein of human ARID2
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.


Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of LO2 cells with ARID2 Polyclonal Antibody
Observed Mw:250kDa
Calculated Mw:191kDa/197kDa



Western blot analysis of extracts from normal (control) and ARID2 knockout (KO) HeLa cells, using ARID2 Polyclonal Antibody at dilution of 1:500.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only