

## ARL8B Polyclonal Antibody

**Catalog No.** E-AB-65602

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                     |                                                          |
|---------------------|----------------------------------------------------------|
| <b>Reactivity</b>   | Human,Mouse                                              |
| <b>Immunogen</b>    | Recombinant fusion protein of human ARL8B (NP_060654.1). |
| <b>Host</b>         | Rabbit                                                   |
| <b>Isotype</b>      | IgG                                                      |
| <b>Purification</b> | Affinity purification                                    |
| <b>Conjugation</b>  | Unconjugated                                             |
| <b>Buffer</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.        |

### Applications Recommended Dilution

|           |            |
|-----------|------------|
| <b>IF</b> | 1:50-1:200 |
|-----------|------------|

### Data



Confocal immunofluorescence analysis of U2OS cells using ARL8B Polyclonal Antibody at dilution of 1:100.  
Blue: DAPI for nuclear staining.

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

ARL8B (ADP Ribosylation Factor Like GTPase 8B) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include GTP binding and obsolete signal transducer activity. An important paralog of this gene is ARL8A.

### For Research Use Only