

ARL8B Polyclonal Antibody

Catalog No. E-AB-65602

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human ARL8B (NP_060654.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF	1:50-1:200
-----------	------------

Data



Confocal immunofluorescence analysis of U2OS cells using ARL8B Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

ARL8B (ADP Ribosylation Factor Like GTPase 8B) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include GTP binding and obsolete signal transducer activity. An important paralog of this gene is ARL8A.

For Research Use Only