

SUPT5H Polyclonal Antibody

Catalog No. E-AB-63943

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant protein of human SUPT5H
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:200-1:2000

Data



Western blot analysis of extracts of various cell lines
with SUPT5H Polyclonal Antibody

Observed Mw:170kDa

Calculated Mw:120kDa/121kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

SUPT5H, also named as DSIF p160, is a 1087 amino acid protein, which belongs to the SPT5 family. SUPT5H localizes in the nucleus and is ubiquitously expressed. SUPT5H is a component of the DRB sensitivity-inducing factor complex (DSIF complex), which regulates mRNA processing and transcription elongation by RNA polymerase II. SUPT5H is modified through phosphorylation and methylated, so the molecular weight of SUPT5H might be 140-150 kDa.

For Research Use Only