

RAB11FIP1 Polyclonal Antibody

Catalog No. E-AB-63944

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat, Monkey
Immunogen	Recombinant protein of human RAB11FIP1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:200 - 1:2000

Data



Western blot analysis of extracts of various cell lines
with RAB11FIP1 Polyclonal Antibody

Observed Mw: 80kDa

Calculated

Mw: 15kDa/47kDa/65kDa/70kDa/137kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene encodes one of the Rab11-family interacting proteins (Rab11-FIPs), which play a role in the Rab-11 mediated recycling of vesicles. The encoded protein may be involved in endocytic sorting, trafficking of proteins including integrin subunits and epidermal growth factor receptor (EGFR), and transport between the recycling endosome and the trans-Golgi network. Alternative splicing results in multiple transcript variants. A pseudogene is described on the X chromosome.

For Research Use Only