

HECTD1 Polyclonal Antibody

Catalog No. E-AB-65624

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant fusion protein of human HECTD1 (NP_056197.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using HECTD1 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:300kDa
Calculated Mw:289kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

HECTD1 (HECT Domain E3 Ubiquitin Protein Ligase 1) is a Protein Coding gene. Diseases associated with HECTD1 include Spastic Paraplegia 32, Autosomal Recessive and Prostate Small Cell Carcinoma. Among its related pathways are Innate Immune System and Class I MHC mediated antigen processing and presentation. Gene Ontology (GO) annotations related to this gene include binding and ubiquitin-protein transferase activity. An important paralog of this gene is HECTD4.

For Research Use Only