

## ZRANB3 Polyclonal Antibody

**Catalog No.** E-AB-65632

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Human, Mouse, Rat
<b>Immunogen</b>	Recombinant fusion protein of human ZRANB3 (NP_001273497.1).
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Applications Recommended Dilution

**WB** 1:1000-1:3000

### Data



Western blot analysis of extracts of various cell lines using ZRANB3 Polyclonal Antibody at dilution of 1:3000.

**Observed Mw:123kDa**

**Calculated Mw:42kDa/49kDa/67kDa/123kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

ZRANB3 (Zinc Finger RANBP2-Type Containing 3) is a Protein Coding gene. Diseases associated with ZRANB3 include Warburg Micro Syndrome and Warburg Micro Syndrome 1. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and endonuclease activity. An important paralog of this gene is CHD7.

### For Research Use Only