

FTSJ3 Polyclonal Antibody

Catalog No. E-AB-62566

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant protein of human FTSJ3
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500-1:1000

Data



Western blot analysis of extracts of various cell lines
with FTSJ3 Polyclonal Antibody

Observed Mw:125kDa

Calculated Mw:96kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli FtsJ protein and to a mouse protein possibly involved in embryogenesis.

For Research Use Only