

ANXA5 Polyclonal Antibody

Catalog No. E-AB-40200

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	Recombinant Human Annexin A5 protein
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western Blot analysis of HeLa, A431, HepG2, k562
and Raw264.7 cells using ANXA5 Polyclonal
Antibody at dilution of 1:500

Observed Mw:36 kDa
Calculated Mw:36 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Annexin A5 (ANXA5), is a member of the structurally related family of annexin proteins some of which have been implicated in membrane-related events along exocytotic and endocytotic pathways. Annexin 5 is a phospholipase A2 and protein kinase C inhibitory protein with calcium channel activity and a potential role in cellular signal transduction, inflammation, growth and differentiation. Annexin 5 has also been described as placental anticoagulant protein I, vascular anticoagulant- α , endonexin II, lipocortin V, placental protein 4 and anchorin CII.

For Research Use Only