

SARS-CoV-2 NSP5 Polyclonal Antibody(2019-nCoV)

Catalog Number:E-AB-V1040



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	SARS-COV2
Immunogen	Recombinant SARS-CoV-2 3C-like Proteinase
Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Formulation	PBS, pH7.4, containing 0.05% proclin300, 50% glycerol.

Applications Recommended Dilution

ELISA	1:4000~1:8000
--------------	---------------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

Cleaves the C-terminus of replicase polyprotein at 11 sites (By similarity). Recognizes substrates containing the core sequence [ILMVF]-Q-I-[SGACN]. Also able to bind an ADP-ribose-1"-phosphate (ADRP) (By similarity).

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017