

# SARS-CoV-2 NSP9 Polyclonal Antibody(2019-nCoV)

Catalog Number:E-AB-V1042



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	SARS-COV2
<b>Immunogen</b>	Recombinant SARS-CoV-2 NSP9
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS, pH7.4, containing 0.05% proclin300, 50% glycerol.

## Applications Recommended Dilution

<b>ELISA</b>	1:4000~1:8000
--------------	---------------

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

## Background

May participate in viral replication by acting as a ssRNA-binding protein.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017