

SARS-CoV-2 NSP12 Polyclonal Antibody(2019-nCoV)

Catalog Number:E-AB-V1044



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	SARS-COV2
Immunogen	Recombinant SARS-CoV-2 NSP12
Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Formulation	PBS, pH7.4, containing 0.05% proclin300, 50% glycerol.

Applications Recommended Dilution

ELISA	1:4000~1:8000
WB	1:500~1000

Data



Western Blot analysis of pseudovirion using SARS-CoV-2 NSP12 polyclonal antibody(2019-nCoV) at dilution of 1:1000.

Observed Mw:100kDa
Calculated Mw:100kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Responsible for replication and transcription of the viral RNA genome.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017