

SARS-CoV-2 ORF8 Protein Polyclonal Antibody (2019-nCoV)



Catalog Number:E-AB-V1048

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	SARS-COV2
Immunogen	Recombinant SARS-CoV-2 ORF8
Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Formulation	PBS, pH7.4, containing 0.05% proclin300, 50% glycerol.

Applications Recommended Dilution

ELISA	1:4000~1:8000
--------------	---------------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

May play a role in host-virus interaction.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017