

## DICER1 Polyclonal Antibody

**Catalog No.** E-AB-70297

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Human,Rat
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Mouse Dicer
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Buffer</b>	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4


### Applications Recommended Dilution

**IHC** 1:300-1:1000

### Data



Immunohistochemistry analysis of paraffin-embedded human kidney using DICER1 Polyclonal Antibody at dilution of 1:400.



Immunohistochemistry analysis of paraffin-embedded rat brain using DICER1 Polyclonal Antibody at dilution of 1:400.

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

This gene encodes a protein possessing an RNA helicase motif containing a DEXH box in its amino terminus and an RNA motif in the carboxy terminus. The encoded protein functions as a ribonuclease and is required by the RNA interference and small temporal RNA (stRNA) pathways to produce the active small RNA component that represses gene expression. Alternative splicing results in multiple transcript variants.

### For Research Use Only