

# Recombinant HCoV-OC43 Nucleocapsid Protein

Catalog Number:PKSV030284



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Synonyms</b>	HCoV-OC43 Nucleocapsid Protein;HCoV-OC43 coronavirus NP Protein;HCoV-OC43 np Protein;HCoV-OC43 novel coronavirus Nucleoprotein Protein
<b>Species</b>	HCoV
<b>Expression Host</b>	E.coli
<b>Sequence</b>	Met1-Ile448
<b>Accession</b>	P33469
<b>Calculated Molecular Weight</b>	53.0 kDa
<b>Observed molecular weight</b>	57 kDa
<b>Tag</b>	N-His

## Properties

<b>Purity</b>	> 90 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	Please contact us for more information.
<b>Storage</b>	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
<b>Shipping</b>	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
<b>Formulation</b>	Supplied as a 0.2 µm filtered solution of 20mM Acetate, 10% Sucrose, 50mM NaCl, 0.05% Tween80, pH4.5
<b>Reconstitution</b>	Not Applicable

## Background

Coronavirus N protein is required for coronavirus RNA synthesis, and has RNA chaperone activity that may be involved in template switch. Nucleocapsid protein is a most abundant protein of coronavirus. N protein packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP) and plays a fundamental role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral replication. Because of the conservation of N protein sequence and its strong immunogenicity, the N protein of coronavirus is chosen as a diagnostic tool.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017