

FLI1 Monoclonal Antibody

Catalog No. E-AB-71008

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length protein from eukaryotic expression system.
Host	Mouse
Isotype	IgG
Clone	LE8A1
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.


Applications Recommended Dilution

IHC 1:100-1:200

Data



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human liver using FLI1 Monoclonal Antibody with Cat#E-AB-71008 at dilution of 1:200.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human tonsil placenta using with Cat#E-AB-71008 at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

FLI1 is a transcription factor containing an ETS DNA-binding domain. It can translocate with the Ewing sarcoma gene on chromosome 22, which results in a fusion gene that is present in the majority of Ewing sarcoma cases. An acute lymphoblastic leukemia-associated translocation involving this gene has also been identified. Diseases associated with FLI1 include Isolated Delta-Storage Pool Disease and Ewing Sarcoma. Among its related pathways are Transcriptional misregulation in cancer and NF-kappaB Signaling.

For Research Use Only