

Oct-2 Monoclonal Antibody

Catalog No. E-AB-71011

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length protein from eukaryotic expression system.
Host	Mouse
Isotype	IgG
Clone	LE0202
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

IHC 1:100-1:200

Data



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human tonsil tissue using Oct-2 Monoclonal Antibody with Cat#E-AB-71011 at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Oct-2 is a transcription factor of the POU family, which can bind to the Ig gene octamer and regulate the expression of the B cell specific gene. Oct-2 protein is expressed in B cells, which is another useful antibody for identifying B cells and B cell lymphomas.

For Research Use Only