

CD31 Monoclonal Antibody

Catalog No. E-AB-71031

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length protein from eukaryotic expression system.
Host	Mouse
Isotype	IgG
Clone	LE0034
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.


Applications Recommended Dilution

IHC 1:100-1:200

Data



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human tonsil tissue using CD31 Monoclonal Antibody with Cat#E-AB-71031 at dilution of 1:200.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human liver tissue using with Cat#E-AB-71031 at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CD31 is mainly applied in marking mononuclear cells, granule cells, endothelial cells and some T cells, and is mostly used for diagnosis and differential diagnosis of benign and malignant vascular tumors, and it can also be used to research the angiogenesis during a variety of tumor interstitial.

For Research Use Only