

CK-13 Monoclonal Antibody

Catalog No. E-AB-71040

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Human
Immunogen	Full length protein from eukaryotic expression system.
Host	Mouse
Isotype	IgG
Clone	LE0075
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Applications Recommended Dilution

IHC 1:100-1:200

Data



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human tonsil tissue using CK-13 Monoclonal Antibody with Cat#E-AB-71040 at dilution of 1:200.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human lung squamous carcinoma using with Cat#E-AB-71040 at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Cytokeratin 13 (CK13) has a molecular weight of 53 KDa, it is mainly marked on non-keratinized squamous epithelium and can be used for diagnosis of squamous cell carcinoma.

For Research Use Only