

PGP 9.5 Monoclonal Antibody

Catalog No. E-AB-71089

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length protein from eukaryotic expression system.
Host	Mouse
Isotype	IgG
Clone	LE0213
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

IHC 1:100-1:200

Data



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human appendix tissue peripheral nerves using PGP 9.5 Monoclonal Antibody with Cat#E-AB-71089 at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Protein gene product 9.5 (PGP 9.5) isolated from bovine brain has been used for immunization. PGP 9.5 is a pan-neuronal marker. The antibody is a useful aid for the classification of neuroendocrine tumors.

For Research Use Only