

pVHL Monoclonal Antibody

Catalog No. E-AB-71097

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length protein from eukaryotic expression system.
Host	Mouse
Isotype	IgG
Clone	LE0229
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

IHC 1:100-1:200

Data



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of Human liver tissue using pVHL Monoclonal Antibody with Cat#E-AB-71097 at dilution of 1:200.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

VHL is a tumor suppressor gene. VHL disease and most patients of sporadic clear cell renal cell carcinoma often showed mutations of the gene. VHL is mainly used for the study of renal cell carcinoma and other malignant tumors.

For Research Use Only