

# Goat Anti-Mouse IgG(H+L)(Elab Fluor<sup>®</sup> 647 conjugated)

Catalog Number:E-AB-1076



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	Mouse
<b>Immunogen</b>	Mouse IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Elab Fluor <sup>®</sup> 647 (Ex/Em:650nm/665nm)
<b>Formulation</b>	PBS with 0.05% sodium azide and 50% glycerol, PH7.4

## Applications Recommended Dilution

<b>Immunofluorescence</b>	1:100-300
---------------------------	-----------

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017