

Goat Anti-Mouse IgG(H+L)(Elab Fluor[®] 647 conjugated)

Catalog Number:E-AB-1076



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

| | |
|---------------------|-----------------------------------------------------|
| Reactivity | Mouse |
| Immunogen | Mouse IgG |
| Host | Goat |
| Purification | Affinity purification |
| Conjugation | Elab Fluor [®] 647 (Ex/Em:650nm/665nm) |
| Formulation | PBS with 0.05% sodium azide and 50% glycerol, PH7.4 |

Applications Recommended Dilution

| | |
|---------------------------|-----------|
| Immunofluorescence | 1:100-300 |
|---------------------------|-----------|

Preparation & Storage

| | |
|----------------|-----------------------------------------------------------------------------------------|
| Storage | Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light. |
|----------------|-----------------------------------------------------------------------------------------|

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017