

P.gingivalis Monoclonal Antibody

Catalog Number:E-AB-48029



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	P.gingivalis
Immunogen	P.gingivalis
Host	Mouse
Isotype	IgG1
Purification	Protein A/G Purified
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol,pH 7.4

Applications Recommended Dilution

IHC	1:100-200
------------	-----------

Data



Immunohistochemistry of P.gingivalis.

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

P. gingivalis invades gingival epithelial cells in high numbers, in which case both bacteria and epithelial cells survive for extended periods of time. High levels of specific antibodies can be detected in patients harboring P. gingivalis. P. gingivalis infection has been linked to Alzheimer's disease and rheumatoid arthritis.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017