

Zebrafish PPP1R3G Polyclonal Antibody

Catalog Number:E-AB-40563



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Zebrafish
Immunogen	Recombinant human PPP1R3G protein expressed by E.coli
Host	Rabbit
Isotype	IgG
Purification	Antigen Affinity Purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol,pH 7.4

Applications Recommended Dilution

WB	1:500-1000
-----------	------------

Data



Western blot of Zebrafish whole lysates with anti-PPP1R3G rabbit polyclonal antibody at dilution of 1:500.

Observed Mw:30kDa
Calculated Mw:30kDa

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

Glycogen-targeting subunit for protein phosphatase 1 (PP1). Involved in the regulation of hepatic glycogenesis in a manner coupled to the fasting-feeding cycle and distinct from other glycogen-targeting subunits.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017