

NCI-H3122 Cells Complete Medium

Cat. No: CM-0721
Size: 125ml×4

General Information

Product Form	Liquid
Concentration	Ready-to-use
Component	RPMI-1640 (PM150110) 10% FBS 1% P/S (PB180120)
Bacterial detection	Negative
Fungal detection	Negative
Mycoplasma detection	Negative
Endotoxin content	≤5
Storage	2-8°C, Shading light
Shipping	Ice bag
Expiration date	3 months

Notes

1. This product is for research use only.
2. This product is sterilized by 0.22 μ m filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution during the culture.
4. It is not suitable for long time storage at room temperature or 4 °C.
5. This product contains serum, penicillin and streptomycin. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.