

Recombinant Human FGF-6 protein(His Tag)

Catalog Number:PKSH034151



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	HBGF-6; HST-2
Species	Human
Expression Host	E.coli
Sequence	Gly 41-Ile 208
Accession	P10767
Calculated Molecular Weight	19.7 kDa
Observed molecular weight	19 kDa
Tag	C-His

Properties

Purity	> 98 % as determined by reducing SDS-PAGE.
Endotoxin	< 0.1 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 20 mM sodium citrate, 0.2 M NaCl, pH 3.5. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

FGF-6 belongs to the fibroblast growth factor (FGF) family. Members of this family possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. FGF6 plays an important role in the regulation of cell proliferation, cell differentiation, angiogenesis and myogenesis. It is also required for normal muscle regeneration. FGF6 gene displayed oncogenic transforming activity when transfected into mammalian cells.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017