

Recombinant Mouse IL-27 EB13 protein(His Tag)

Catalog Number:PKSM041472



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Lymphopoietin 1(LP-1);pre-B-cell factor
Species	Mouse
Expression Host	E.coli
Sequence	Ala 23-Pro 228
Accession	O35228
Calculated Molecular Weight	23.8 kDa
Observed molecular weight	25 kDa
Tag	C-His
Bioactivity	Measure by its ability to protect HepG2 cells infected with encephalomyocarditis (EMC) virus.The ED ₅₀ for this effect is < 5 ng/mL.

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 0.1 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Possesses antiviral functions and induces the type I interferon response during influenza infection. Signals through two receptor complexes IL20RA/IL20RB or IL20RB/IL22RA1. In turn, stimulates the JAK1-STAT3 and MAPK pathways and promotes the secretion of pro-inflammatory mediators including IL8 and MMP1 (By similarity).

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017