

Recombinant Swine IL-6 protein(His Tag)

Catalog No. PKSS000005

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	IFN- β 2;B-Cell Differentiation Factor (BCDF);BSF-2;HPGF;HSF;MGI-2
Species	Porcine
Expression Host	E.coli
Sequence	Arg 31-Met 212
Accession	P26893
Calculated Molecular Weight	21.9 kDa
Observed molecular weight	23 kDa
Tag	C-His
Bioactivity	Measure by its ability to induce proliferation in T1165.85.2.1 cells.The ED ₅₀ for this effect is < 1.5 ng/mL.

Properties

Purity	> 98 % as determined by reducing SDS-PAGE.
Endotoxin	< 0.1 EU per μ g of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 98 % as determined by reducing SDS-PAGE.

Background

Interleukin-6 (IL-6) is a multifunctional α -helical cytokine that regulates cell growth and differentiation of various tissues,

For Research Use Only

which is known particularly for its role in the immune response and acute phase reactions. IL-6 protein is secreted by a variety of cell types including T cells and macrophages as phosphorylated and variably glycosylated molecule. It exerts actions through the its heterodimeric receptor composed of IL-6R that lacks the tyrosine/kinase domain and binds IL-6 with low affinity, and ubiquitously expressed glycoprotein 130 (gp130) that binds the IL-6. IL-6R complex with high affinity and thus transduces signals. IL-6 is also involved in hematopoiesis, bone metabolism, and cancer progression, and has been defined an essential role in directing transition from innate to acquired immunity.