

Recombinant Human CXCL13/BCA-1 Protein

Catalog Number:PDEH100076



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	C-X-C motif chemokine 13;Angie;B cell-attracting chemokine 1;B lymphocyte chemoattractant;CXC chemokine BLC;Small-inducible cytokine B13;BCA1;BLC;SCYB13
Species	Human
Expression Host	E.coli
Sequence	Val23-Arg94
Accession	O43927
Calculated Molecular Weight	8.7 kDa
Observed molecular weight	10 kDa
Tag	N-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4., 5% trehalose, 5% mannitol, 0.01% tween-80. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the print
Reconstitution	Please refer to the printed manual for detailed information.

Background

Chemotactic for B-lymphocytes but not for T-lymphocytes, monocytes and neutrophils. Does not induce calcium release in B-lymphocytes. Binds to BLR1/CXCR5.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017