

Recombinant Human AZGP1 protein (His tag)

Catalog Number:PDEH100173



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	Zinc-Alpha-2-Glycoprotein;Zn-Alpha-2-GP;Zn-Alpha-2-Glycoprotein;AZGP1;ZAG;ZNGP1
Species	Human
Expression Host	E.coli
Sequence	Asn 148-Ser 298
Accession	P25311
Calculated Molecular Weight	16.5 kDa
Observed molecular weight	22 kDa
Tag	N-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Zinc- α -2-Glycoprotein (AZGP1) can be found in blood plasma; seminal plasma; urine; sweat; saliva; liver; and epithelial cells of various human glands. AZGP1 has been proposed in the regulation of body weight; and the melanin production by normal and malignant melanocytes. AZGP1 stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. AZGP1 has been reported to stimulate lipid breakdown and may have an important role in lipid homeostasis. Mature human AZGP1 consists of one MHC class I antigen region and a C2-type Ig-like domain. AZGP1 has two alternate splice forms; one shows a 66 amino acids substitution for the C-terminal 30 amino acids; the other one shows a nine Lys substitution for amino acid 151-298.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017