

Recombinant Human Interferon α -1/IFNA1 Protein (His Tag)



Catalog Number: PDMH100009

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Interferon alpha-1/13; IFN-alpha-1/13; Interferon alpha-D; LeIF D; IFNA1; IFNA13
Species	Human
Expression Host	HEK293 Cells
Sequence	Cys24-Glu189
Accession	P01562
Calculated Molecular Weight	20.4 kDa
Observed molecular weight	17 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μ g of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Interferon alpha-1/13 (IFN-alpha-1/13 for short), also known as Interferon alpha-D, is a secreted protein which belongs to the alpha/beta interferon family. It is produced by macrophages. IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. IFN-alpha exerts a variety of other biological effects, including antitumor and immunomodulatory activities and are increasingly used clinically to treat a range of malignancies, myelodysplasias and autoimmune diseases.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017