Recombinant Mouse S100A15/S100A7A Protein

Catalog Number:PDEM100003



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description	
Synonyms	S100 calcium-binding protein A15A;Protein S100-A15A;Protein S100-A7A;S100 calcium-binding protein A7A;S100a15a
Species	Mouse
Expression Host	E.coli
Sequence	Met1-Tyr108
Accession	Q6S5I3
Calculated Molecular Weight	12.9 kDa
Observed molecular weight	12 kDa
Tag	None
Properties	
Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μ g of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Koebnerisin is also known as protein S100-A7A (S100A7A), S100 calcium-binding protein A7-like 1 (S100A7L1) or S100 calcium-binding protein A15 (S100A15). Human S100A7A / S100A15 is a novel member of the S100 family of EF-hand calcium-binding proteins and was recently identified in psoriasis, where it is significantly upregulated in lesional skin. S100A7 is expressed by both normal cultured and malignant keratinocytes and malignant breast epithelial cells within ductal carcinoma in situ, suggesting an association with abnormal pathways of differentiation. S100A7 plays a role in the pathogenesis of inflammatory skin disease, as a chemotactic factor for hematopoietic cells. It also plays a role in early stages of breast tumor progression in association with the development of the invasive phenotype. The association of the 11.2 kDa S100A7A / S100A15 with psoriasis suggests that it contributes to the pathogenesis of the disease and could provide a molecular target for therapy.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine Toll-free: 1-888-852-8623 Tel: 1-832-243-608