

ADGRL4 Polyclonal Antibody

Catalog No. E-AB-90900

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human ADGRL4
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using ADGRL4 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:78kDa
Calculated Mw:77kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze/thaw cycles.

Background

Endothelial orphan receptor that acts as a key regulator of angiogenesis.

For Research Use Only