

NUDCD1 Polyclonal Antibody

Catalog No. E-AB-91352

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human NUDCD1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using NUDCD1 Polyclonal Antibody at 1000 dilution.

Observed Mw: 73kDa

Calculated Mw: 56kDa/63kDa/66kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze/thaw cycles.

Background

Predicted to be involved in immune system process. Located in cytosol and nucleoplasm.

For Research Use Only