

NUBP1 Polyclonal Antibody

Catalog No. E-AB-91520

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse
Immunogen	Recombinant fusion protein of human NUBP1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using NUBP1 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:35kDa

Calculated Mw:33kDa/34kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze/thaw cycles.

Background

NUBP1 is a member of the NUBP/MRP subfamily of ATP-binding proteins.

For Research Use Only