

CNTD1 Polyclonal Antibody

Catalog No. E-AB-91920

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human
Immunogen	Recombinant fusion protein of human CNTD1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Applications

Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using CNTD1 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:37kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze/thaw cycles.

Background

Predicted to be involved in reciprocal meiotic recombination. Predicted to act upstream of or within spermatogenesis. Predicted to be active in site of double-strand break.

For Research Use Only