

# Rabbit IgG

Catalog Number:E-AB-1128



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Antigen Affinity Purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	0.02M PBS,pH7.2,50% glycerol

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze/thaw cycles.
----------------	---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017