

Elab Fluor[®] 405-Goat Anti-Rabbit IgG F(ab')₂

Catalog Number:E-AB-1113



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rabbit
Immunogen	Rabbit IgG F(ab') ₂
Host	Goat
Purification	Affinity purification
Conjugation	Elab Fluor [®] 405 (Ex/Em:407nm/488nm)
Formulation	PBS with 0.02% sodium azide,1% protective protein and 50% glycerol,pH7.4

Applications Recommended Dilution

IF	1:150-1:300
FCM	1:250-1:500

Preparation & Storage

Storage	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
Shipping	The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017