

# Biotin-Goat Anti-Rabbit IgG F(ab')<sub>2</sub>

Catalog Number:E-AB-1114



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	Rabbit
<b>Immunogen</b>	Rabbit IgG F(ab') <sub>2</sub>
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Biotin
<b>Formulation</b>	PBS with 0.02% sodium azide and 1% protective protein,pH7.4

## Applications Recommended Dilution

<b>WB</b>	1:5000-1:10000
<b>IHC</b>	1:250-1:500
<b>ELISA</b>	1:5000-1:10000

## Preparation & Storage

<b>Storage</b>	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
<b>Shipping</b>	The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017