

# Elab Fluor<sup>®</sup> 594-Rabbit Anti-Rat IgG(H+L)

Catalog Number:E-AB-1152



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

|                     |  |
|---------------------|--|
| <b>Reactivity</b>   | Rat  |
| <b>Immunogen</b>    | Rat IgG  |
| <b>Host</b>         | Rabbit   |
| <b>Purification</b> | Affinity purification  |
| <b>Conjugation</b>  | Elab Fluor <sup>®</sup> 594 (Ex/Em:591nm/614nm)                          |
| <b>Formulation</b>  | PBS with 0.02% sodium azide,1% protective protein and 50% glycerol,pH7.4 |

## Applications Recommended Dilution

|            |             |
|------------|-------------|
| <b>IF</b>  | 1:150-1:300 |
| <b>FCM</b> | 1:250-1:500 |

## Preparation & Storage

|                 |   |
|-----------------|---|
| <b>Storage</b>  | Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light. |
| <b>Shipping</b> | The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended.      |

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017