

FITC-Rabbit Anti-Mouse IgG

Catalog Number:E-AB-1158



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Rabbit
Purification	Affinity purification
Conjugation	FITC (Ex/Em:495nm/525nm)
Formulation	PBS with 0.02% sodium azide and 1% gelatin,pH7.4

Applications Recommended Dilution

FCM	1:250-1:500
------------	-------------

Preparation & Storage

Storage	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
Shipping	The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017