

Elab Fluor[®] 647-Mouse Anti-Rabbit IgG(H+L)

Catalog Number:E-AB-1180



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rabbit
Immunogen	Rabbit IgG
Host	Mouse
Purification	Affinity purification
Conjugation	Elab Fluor [®] 647 (Ex/Em:650nm/665nm)
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

FCM	1:250-1:500
------------	-------------

Preparation & Storage

Storage	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
Shipping	The product is shipped with ice pack, upon receipt, store it immediately at the temperature recommended.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017