

# Elab Fluor<sup>®</sup> 594-Goat Anti-Guinea Pig IgG (H+L)

Catalog Number:E-AB-1187



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	Guinea Pig
<b>Immunogen</b>	Guinea Pig IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Elab Fluor <sup>®</sup> 594 (Ex/Em:591nm/614nm)
<b>Formulation</b>	PBS with 0.02% sodium azide,1% protective protein and 50% glycerol,pH7.4

## Applications Recommended Dilution

<b>IF</b>	1:150-1:300
<b>FCM</b>	1:250-1:500

## Preparation & Storage

<b>Storage</b>	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
<b>Shipping</b>	The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017