

Elab Fluor[®] 405-Goat Anti-Guinea Pig IgG (H+L)

Catalog Number:E-AB-1189



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

| | |
|---------------------|--|
| Reactivity | Guinea Pig |
| Immunogen | Guinea Pig IgG |
| Host | Goat |
| Purification | Affinity purification |
| Conjugation | Elab Fluor [®] 405 (Ex/Em:407nm/488nm) |
| Formulation | PBS with 0.02% sodium azide,1% protective protein and 50% glycerol,pH7.4 |

Applications Recommended Dilution

| | |
|------------|-------------|
| IF | 1:150-1:300 |
| FCM | 1:250-1:500 |

Preparation & Storage

| | |
|-----------------|---|
| Storage | Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light. |
| Shipping | The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended. |

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017