

HRP-Donkey Anti-rabbit IgG(H+L) preadsorbed

Catalog Number:E-AB-1193



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rabbit
Immunogen	rabbit IgG
Host	Donkey
Purification	Affinity purification
Conjugation	HRP
Formulation	PBS with 0.1% procline 300,1% protective protein and 50% glycerol,pH7.4

Applications Recommended Dilution

WB	1:5000-1:10000
IHC	1:250-1:500
ELISA	1:5000-1:10000

Preparation & Storage

Storage	Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
Shipping	The product is shipped with ice pack,upon receipt,store it immediately at the temperature recommended.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017